


# Manual

KGO Enzyme Panel

**AminoAcid**

Fe(II)/ $\alpha$ -ketoglutarate-dep. oxygenases



Enzyme panel of KGOs for  
**Oxyfunctionalization**

# Content

1.	Introduction	3
2.	KGO Enzyme Panel	4
3.	AminoAcid – KGO Enzyme Panel	5
4.	Example reactions	6
5.	Components and storage	7
6.	Plate content	8
7.	Plate layout	9
8.	Required reagents and equipment	10
9.	Preparation of reagents	11
10.	Protocol	13
11.	Pipetting scheme	14
12.	Adapting the protocol	15
13.	Miscellaneous	16


## 1

# Introduction

Fe(II)- and  $\alpha$ -ketoglutarate dependent oxygenases (KGOs) are non-heme iron proteins, that rely on  $\alpha$ -ketoglutarate as a cofactor. These enzymes catalyse a variety of interesting reactions including, among others, hydroxylation, epoxidation, demethylation, ring formation or ring expansion.

---

**$\alpha$ -KGOs are ideal biocatalysts for chemically challenging reactions**




## 2

## KGO Enzyme Panel

The KGO Enzyme Panel contains wild-type KGOs recombinantly produced in *Escherichia coli*.

The practical MTP format allows parallel screening of the entire KGO Enzyme Panel, followed by your GC, HPLC or MS analysis of choice.

Any KGO from the KGO Enzyme Panel yielding the desired product can be engineered towards even greater performance.



## 3

## AminoAcid – KGO Enzyme Panel

The AminoAcid KGO Enzyme Panel contains **computationally curated** wild-type KGOs recombinantly produced in *Escherichia coli* that act on amino acids such as the following:

- L-Leucine
- L-Proline
- ...

Why curate?

Due to the inherently **higher selectivity of KGOs** compared to UPOs for example, the curation reduces overall panel size and thus screening efforts while simultaneously **increasing chances to identify a hit**.

Aminoverse also offers 3 other curated KGO Enzyme Panels focused on the modification of different metabolites:


1. **BioCarbon**      fatty acids, terpenes, steroids
2. **XNA**              nucleic acids, nucleosides and –tides
3. **Eureka**            discovery panel with KGOs not fitting into the other three panels

Contact us for more information!  
[enzymes@aminoverse.com](mailto:enzymes@aminoverse.com)

## 4

## Example reactions

Examples for amino acid hydroxylation

[1] P. Lukat et al., *Chem. Sci.* **2017**, 8, 7521–7527[2] T. J. Smart et al., *Bioorg. Chem.* **2020**, 94[3] X. Yin, T. Zabriskie, *ChemBioChem* **2004**, 5, 1274–1277

## 5

## Components and storage

### KGO Enzyme Panel

- 1x MTP, contains lyophilized *E. coli* cell lysate supernatant + negative control (see plate layout on p. 9)
- Store at -20 °C
- Use directly after resuspension
- Caution: avoid repeated freeze-thaw cycles of the KGO Enzyme Panel



The KGO Enzyme Panel MTP contains:

- 42 KGOs in duplicates
- Negative control in duplicates (well C6, C12)
- Extra KGO 11 in duplicates (well D6, D12), which can be used as positive control for test reactions to get familiar with the handling and the reaction before the actual screening

Reactions possible with 1 KGO Enzyme Panel:

- If you follow the recommended reaction protocol you can perform 1x 200  $\mu$ L reactions per KGO



7

Plate layout

1 2 3 4 5 6 7 8 9 10 11 12

A	KG0 28	KG0 49	KG0 75	KG0 88	KG0 100	KG0 246	KG0 28	KG0 49	KG0 75	KG0 88	KG0 100	KG0 246
B	KG0 30	KG0 50	KG0 77	KG0 90	KG0 102	KG0 247	KG0 30	KG0 50	KG0 77	KG0 90	KG0 102	KG0 247
C	KG0 32	KG0 57	KG0 78	KG0 93	KG0 131	neg. c. (no KGO)	KG0 32	KG0 57	KG0 78	KG0 93	KG0 131	neg. c. (no KGO)
D	KG0 34	KG0 60	KG0 80	KG0 94	KG0 146	KG0 11	KG0 34	KG0 60	KG0 80	KG0 94	KG0 146	KG0 11
E	KG0 40	KG0 61	KG0 81	KG0 95	KG0 217		KG0 40	KG0 61	KG0 81	KG0 95	KG0 217	
F	KG0 41	KG0 66	KG0 84	KG0 96	KG0 218		KG0 41	KG0 66	KG0 84	KG0 96	KG0 218	
G	KG0 45	KG0 69	KG0 85	KG0 98	KG0 221		KG0 45	KG0 69	KG0 85	KG0 98	KG0 221	
H	KG0 47	KG0 74	KG0 87	KG0 99	KG0 222		KG0 47	KG0 74	KG0 87	KG0 99	KG0 222	

## 8

## Required reagents and equipment

- (Multi-channel) micro pipette
- dH<sub>2</sub>O
- Buffer (*e.g.* 1000 mM KP<sub>i</sub> buffer, pH 6.3)
- Sodium ketoglutarate
- Sodium ascorbate
- FeSO<sub>4</sub>
- Substrate stock solution
- Reaction vessels, such as *e.g.*:
  - Microtiter plates
  - Reaction tubes
  - Glass vials
- Shaking incubator for 25 °C and 500 rpm

## 9

## Preparation of reagents

- **KGO Enzyme Preparation**
  - Resuspend lyophilized KGO Enzyme Preparation/negative control with 115  $\mu\text{L}$   $\text{dH}_2\text{O}$  per well
- **Buffer**
  - Prepare buffer of choice, e.g. 1000 mM  $\text{KP}_i$  buffer, pH 6.3
- **Sodium ketoglutarate stock**
  - Prepare 100 mM stock in  $\text{dH}_2\text{O}$
  - 0.9 mL\* are sufficient for testing 42 KGOs plus negative control in 200  $\mu\text{L}$  reactions
- **Sodium ascorbate stock**
  - Prepare 5 mM stock in  $\text{dH}_2\text{O}$
  - 1.2 mL\* are sufficient for testing 42 KGOs plus negative control in 200  $\mu\text{L}$  reactions

Turn page for the rest of preparations



\* Sufficient volume to accommodate for pipetting errors

- **FeSO<sub>4</sub> stock**

- Prepare 50 mM stock in dH<sub>2</sub>O
- 1.2 mL\* are sufficient for testing 42 KGOs plus negative control in 200 µL reactions

- **Substrate stock**

- Prepare 50 mM substrate stock solution, *e.g.* dissolved in dH<sub>2</sub>O
- 1.2 mL\* are sufficient for testing 42 KGOs plus negative control in 200 µL reactions

**Note:** Instead of adding the components individually, a master mix, containing KP<sub>i</sub> buffer, sodium ketoglutarate, sodium ascorbate, FeSO<sub>4</sub> and substrate can be prepared and added to the resuspended KGOs accordingly.

\* Sufficient volume to accommodate for pipetting errors


## 10 Protocol

1. Centrifuge plate before opening it
2. Add **115  $\mu\text{L}$  of resuspended KGO solution** or resuspended negative control to suitable reaction vessel, *e.g.* microtiter plate
3. Add **10  $\mu\text{L}$  of 1000 mM  $\text{KPi}$  buffer pH 6.3** to reaction vessel
4. Add **15  $\mu\text{L}$  of 100 mM sodium ketoglutarate** stock solution
5. Add **20  $\mu\text{L}$  of 5 mM sodium ascorbate** stock solution
6. Add **20  $\mu\text{L}$  of 50 mM  $\text{FeSO}_4$**  stock solution
7. Add **20  $\mu\text{L}$  of 50 mM substrate** stock solution to start the reaction
7. Incubate at **25  $^\circ\text{C}$  for 4 h\* at 500 rpm**
8. Continue with your preferred extraction protocol, *e.g.* ethyl acetate extraction and drying over  $\text{Na}_2\text{SO}_4$  to perform GC analysis

\* Extending reaction time up to 24 h can lead to higher conversion rates

Compound	Volume [μL]	Stock conc.	Final conc.
Resuspended KGO solution/ negative control	115	–	–
KP <sub>i</sub> Buffer, pH 6.3*	10	1000 mM	50 mM
Sodium ketoglutarate*	15	100 mM	7.5 mM
Sodium ascorbate*	20	5 mM	0.5 mM
FeSO <sub>4</sub> × 7 H <sub>2</sub> O*	20	50 mM	5 mM
Substrate solution*	20	50 mM	5 mM
Total volume	200		

- Reactions can be performed e.g. in microtiter plates or reaction tubes



\* **Note:** Instead of adding the components individually, a master mix, containing KP<sub>i</sub> buffer, sodium ketoglutarate, sodium ascorbate, FeSO<sub>4</sub> and substrate can be prepared and added to the resuspended KGOs accordingly.

- Reaction volumes, reagent concentrations, and incubation times may be adapted to the desired reaction and the sensitivity of the detection method
- Reaction solutions can be analyzed using GC, HPLC or MS. Depending on the chosen analysis method, the extraction procedure may vary

### Reaction volume

- Can be adjusted to specific needs; *e.g.* scaled up to 500  $\mu\text{L}$

### Reaction time

- Can be varied according to specific needs and conversion of the substrate

### Substrate stock solution and cosolvents

- Substrate concentration can be adjusted to specific needs and according to solubility

Feel free to contact us in case of questions or for tuning of the reaction conditions at [enzymes@aminoverse.com](mailto:enzymes@aminoverse.com).

The KGO Enzyme Panel is for *in vitro* research use only, not for diagnostic or therapeutic applications and has to be handled by qualified personnel.

The purchase of this product only authorizes the buyer to application of the KGO enzyme panel for internal research. Consult Aminoverse's Terms and Conditions for more information.

Reference information:

In publications, please refer to the KGOs of this enzyme panel by their *number* as “Aminoverse KGOxxx”.





Aminoverse solves enzyme challenges.

Founded in 2020, Aminoverse offers innovative enzyme products and services:

- Enzyme products and kits  
Ready-to-use enzymes for proof-of-concept studies up to commercial scale, e.g. the KGO Enzyme Panel or the UPO Enzyme Panel, and analyte detection kits, e.g. the Phosfinity series.
- CRO services  
Discovery of enzymes, enzyme feasibility studies, engineering of enzymes by Directed Evolution and machine learning, enzyme mutant library generation.

You have questions regarding this manual or about the application of the KGO enzyme panel?

We look forward to hearing from you.

**Aminoverse B.V.**

Daelderweg 9

6361HK Nuth

The Netherlands

[enzymes@aminoverse.com](mailto:enzymes@aminoverse.com)

[www.aminoverse.com](http://www.aminoverse.com)

+31 4520 848 19

**Aminoverse B.V.**

Daelderweg 9  
6361HK Nuth  
The Netherlands

[enzymes@aminoverse.com](mailto:enzymes@aminoverse.com)

[www.aminoverse.com](http://www.aminoverse.com)

+31 4520 848 19

